



Custom Seating & Positioning Products

# Torso Belt Order Form

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Resale No: \_\_\_\_\_

Contact: \_\_\_\_\_ Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Bill To: \_\_\_\_\_ Ship To: \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_

**\*\* Each Order is Required to Provide Measurements for Users Chest Size, Chair Circumference & Chair Width\*\***

## 1. What Size Chest Strap?

- A) For Individual Users; Give Chest Circumference in Inches: \_\_\_\_\_
- B) Designate Size: Small (28-33") Medium (30-38") Large (36-47") Extra Large (44-53")

## 2. What Material Will be Used for the Chest Strap?

- A) 4" Loop - Most Commonly Ordered
- B) 2" Webbing - Recommended Only if Chaffing is Not a Concern

## 3. How Will the Torso Support Open/Close?

- A) Closure 1 - Hook and Loop Overlap, With D-Ring
- B) Closure 2 - Thread Through and Fold Back, With D-Ring
- C) Closure 3 - Additional \$9.00 for Plastic Quick - Release Buckle, or Additional \$ 14.00 for Auto-Style Buckle, or \$18.00 for Airline Buckle

Note: Multi Adjustable for Closures 1 & 2

## 4. How Will the Webbing Back Strap Attach?

### A) Buckle in Back

Circumference of Back: \_\_\_\_\_ (Always Needed)

Width of Back: \_\_\_\_\_ (Always Needed)

### B) Attach with Grommets - Does Not Apply with Recline Systems

Width of Chair: \_\_\_\_\_

Distance from Screw to Screw; Across the Front- Following the Contour - Including Back Cushion or Laterals: \_\_\_\_\_

6" Torso Belt Extensions are Also Available, (Applicable for Closure 2 Belts) Circle if Wanted -- Retail \$25.96

**Note: The Torso Belt is not compatible with all recline systems. There must be a location to attach the Webbing Back Strap that will allow the Torso Belt to move with the User when reclining.**

**\*\*Torso Belts Retail at \$250 ; Dealer Discounts Vary \*\***

***This product is a Positioning Belt only. It is intended to be used as a mild restraint. It is not designed for use as a safety device that requires high stress loads. If signs of wear appear, the Torso Belt must be replaced.***

Nutec Rehab Inc.  
5840 S. Watt Ave. Suite A  
Sacramento, CA 95829

P 888-448-0093  
F 833-722-0013  
EM customerservice@nutecseating.com