



Custom Seating & Positioning Products

Torso Belt Order Form

Date: _____ PO#: _____ Resale No: _____

Contact: _____ Phone: _____ Fax: _____

Bill To: _____ Ship To: _____

**** Each Order is Required to Provide Measurements for Users Chest Size, Chair Circumference & Chair Width****

1. What Size Chest Strap?

- A) For Individual Users; Give Chest Circumference in Inches: _____
- B) Designate Size: Small (28-33") Medium (30-38") Large (36-47") Extra Large (44-53")

2. What Material Will be Used for the Chest Strap?

- A) 4" Loop - Most Commonly Ordered
- B) 2" Webbing - Recommended Only if Chaffing is Not a Concern

3. How Will the Torso Support Open/Close?

- A) Closure 1 - Hook and Loop Overlap, With D-Ring
- B) Closure 2 - Thread Through and Fold Back, With D-Ring
- C) Closure 3 - Additional \$9.00 for Plastic Quick - Release Buckle, or Additional \$ 14.00 for Auto-Style Buckle, or \$18.00 for Airline Buckle

Note: Multi Adjustable for Closures 1 & 2

4. How Will the Webbing Back Strap Attach?

- A) Buckle in Back
 - Circumference of Back: _____ (Always Needed)
 - Width of Back: _____ (Always Needed)
- B) Attach with Grommets - Does Not Apply with Recline Systems
 - Width of Chair: _____
 - Distance from Screw to Screw; Across the Front- Following the Contour - Including Back Cushion or Laterals: _____

6" Torso Belt Extensions are Also Available, (Applicable for Closure 2 Belts) Circle if Wanted -- Retail \$30.00

Note: The Torso Belt is not compatible with all recline systems. There must be a location to attach the Webbing Back Strap which allows the Torso Belt to move with the User when reclining.

****Torso Belts Retail at \$205; Dealer Discounts Vary****

This product is a Positioning Belt only. It is intended to be used as a mild restraint. It is not designed for use as a safety device that requires high stress loads. If signs of wear appear, the Torso Belt must be replaced.